


Product Name: Recombinant Human SRP68 Protein
Catalog #: PEH9888



Summary

Name	SRP68_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q9UHB9
Host	E.coli
Species	Human
Molecular Mass	59 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GGGGGSGGGG GSGGGGSGGG RGAGGEENKE NERPSAGSKA NKEFGDSLSEILQIIKESQ QQHGLRHGDF QRYRGYCSRR QRRLRKTLEF KMGNRHKFTGKKVTEELLTD NRYLLLVLMQ AERAWSYAMQ LKQEANTEPR

image



Background

Alternative Names	Signal recognition particle subunit SRP68 (SRP68) (Signal recognition particle 68 kDa protein)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use