


Product Name: Recombinant Human SNX29 Protein
Catalog #: PEH9892



Summary

Name	SNX29_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TEQ0
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RQATVAMMNR KDELEENRS LRNLLDGEME HSAALRQEVD TLKRKVAEQE ERQGMKVQAL ARENEVLKVQ LKKYVGAVQM LKREGQTAEV PNLWSVDGEV TVAEQKPGEI AEELASSYER KLIEVAEMHG ELIEFNERLH RALVAKEALV SQMRQELIDL RGPVPGDLSQ TSEDQSLSDF

image



Background

Alternative Names Sorting nexin-29 (RUN domain-containing protein 2A)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human SNX29 Protein
Catalog #: PEH9892

