



Product Name: Recombinant Human SPATC1 Protein
Catalog #: PEH9893



Summary

Name	SPERI_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q76KD6
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FNTSDTQAQP SAAQEQVPA SVPTSPTTSP TVTVLASAPA LAPQVATSYT PSSTTHIAQG APHPPSRMHN SPTQNLVPVPH CPPHNAHSPP RTSSSPASVN DSRGPRTTEP STKSMMMEVER KLAHRKTSKF PENPRESKQL

image



Background

Alternative Names	Speriolin (Spermatogenesis and centriole-associated protein 1) (Spermatogenesis-associated protein 15) (Spermatogenic cell-specific Cdc20-binding protein)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human SPATC1 Protein
Catalog #: PEH9893

