


Product Name: Recombinant Human SPATC1L Protein
Catalog #: PEH9894



Summary

Name	SPC1L_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H0A9
Host	E.coli
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QVRLLKENQM LRRLLSQSCQ EGGGHDLLPP RAHAYPEAGS PGSGVPDFGR FTSVADTPSQ LQTSSLEDLL CSHAPLSSD DTSPGCAAPS QAPFKAFLSP PEPHSHRGTD RKLSPLLSPL QDSLVDKTL EPREMVRPKK VCFSESSLPT GDRTRRSYYL NEIQSFAGAE KDARVVEIA FQLDRRILAY VFPGVTRLYG FTVANIPEKI EQTSTKSLDG SVDERKLREL TQRYLALSAR LEKLGYSRDV HPAFSEFLIN TYGILKQRPD LRANPLHSSP AALRKLVIDV VPPKFLGDSL LLLNCLCELS

image



Background

Alternative Names	Speriolin-like protein (Spermatogenesis and centriole-associated protein 1-like protein)
--------------------------	--

Product Name: Recombinant Human SPATC1L Protein
Catalog #: PEH9894



Note

For Research Use Only , Not for Diagnostic Use