


Product Name: Recombinant Human SPAG5 Protein
Catalog #: PEH9896



Summary

Name	SPAG5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96R06
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RSRQLTEKLT VKSQQALQER DVAIEEKQEV SRVLEQVSAQ LEECKGQTEQ LELENSRLAT DLRAQLQILA NMDSQLKELQ SQHTHCAQDL AMKDELLCQL TQSNEEQAAQ WQKEEMALKH MQAELQQQQA VLAKEVRDLK ETLEFADQEN QVAHLELGQV ECQLKTTLEV LRERSLQCEN LKDTVENLTA KLASTIADNQ EQDLEKTRQY SQKLGLLTEQ LQSLTLFLQT KLKEKTEQET LLLSTACPPT QEHPLPNDRT

image



Background

Alternative Names Sperm-associated antigen 5 (Astrin) (Deepest) (Mitotic spindle-associated protein p126) (MAP126)

Product Name: Recombinant Human SPAG5 Protein
Catalog #: PEH9896



Note

For Research Use Only , Not for Diagnostic Use