


Product Name: Recombinant Human SPATA19 Protein
Catalog #: PEH9898



Summary

Name	SPT19_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q7Z5L4
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SSDIDVVESE AVSVLHHWLK KTEEEASRGI KEKLSINHPS QGVREKMSTD SPPTHGQDIH VTRDVVKHHL SKSDLLANQS QEVLEERTRI QFIRWSHTRI FQVPSEMTED IMRDRIEQVR RSISRLTDVS AQDFSMRPSS SDC

image



Background

Alternative Names	Spermatogenesis-associated protein 19, mitochondrial (Spermatogenic cell-specific gene 1 protein) (Spergen-1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use