


Product Name: Recombinant Human VIPAS39 Protein
Catalog #: PEH9901



Summary

Name	SPE39_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H9C1
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YWNSSKFKAF TFDDEDELS QLKESKRAVN SLRDFVDDDD DDDLERSWS GEPVGSISWS IRETAGNSGS THEGREQLKS RNSFSSYAQL PKPTSTYSLS SFFRGRTRPG SFQSLSDALS DTPAKSYAPE LGRPKEGYRD YSNDWSPSDT VRRLRKGVK SLERFRSLQD KLQLLEEAVS

image



Background

Alternative Names	Spermatogenesis-defective protein 39 homolog (hSPE-39) (VPS33B-interacting protein in apical-basolateral polarity regulator) (VPS33B-interacting protein in polarity and apical restriction)
--------------------------	--

Note

Product Name: Recombinant Human VIPAS39 Protein
Catalog #: PEH9901



For Research Use Only , Not for Diagnostic Use