


Product Name: Recombinant Human SPRYD7 Protein
Catalog #: PEH9902



Summary

Name	SPRY7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5W111
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	CCRDGGTGHI PLKEMPAVQL DTQHMGTDVV IVKNGRRICG TGGCLASAPL HQNKSYPEFK IQSTGIWGIG VATQKVNLNQ IPLGRDMHSL VMRNDGALYH NNEEKNRLPA NSLPQEGDVV GITYDHVELN VYLNGKNMHC PASGIRGTVY PVVYVDDSAI LDCQFSEFYH TPPPGEKIL

image



Background

Alternative Names	SPRY domain-containing protein 7 (Chronic lymphocytic leukemia deletion region gene 6 protein) (CLL deletion region gene 6 protein)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human SPRYD7 Protein
Catalog #: PEH9902

