

Summary

Name	SUCHY_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9HAC7
Host	E.coli
Species	Human
Molecular Mass	45 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QSDMNNIKPL EGVKILDTR VLAGPFATMN LGDLGAEVIK VERPGAGDDT RTWGPFPVGT ESTYYLSVNR NKKSIAVNIK DPKGVKIIKE LAAVCDVFVE NYVPGKLSAM GLGYEDIDEI APHIIYCSIT GYGQTGPISQ RAGYDAVASA VSGLMHITGP ENGDVVRPGV AMTDLATGLY AYGAIMAGLI QKYKTGKGLF IDCNLLSSQV ACLSHIAANY LIGQKEAKRW GTAHGSIVPY QAFKTKDGYI VVGAGNNQQF ATVCKILDLP ELIDNSKYKT NHLRVHNRKE LKILSERFE EELTSKWLYL FEGSGVPYGP INNMKNVFAE PQVLHNGLVM EMEHPTVGKI SVPGPAVRYS KFKMSEARPP PLLGQHTTHI LKEVLRYYDDR AIGELLSAGV VDQHETH

image



Background

Product Name: Recombinant Human SUGCT Protein
Catalog #: PEH9908



Alternative Names

Succinate--hydroxymethylglutarate CoA-transferase (EC 2.8.3.13) (Dermal papilla-derived protein 13) (SuccinylCoA:glutarate-CoA transferase)

Note

For Research Use Only , Not for Diagnostic Use