


Product Name: Recombinant Human SAPCD1 Protein
Catalog #: PEH9910



Summary

Name	SAPC1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5SSQ6
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LVQAPYTVLL LPLGTSRQDP GAQSFFLWLR RMQALEREQD ALWQGLELLQ HGQAWFEDHL REAQRQQQLHL GALGENFLTD LHSEPGRPPL AQIQKVNICL QNLIHEKELS RQQKGVTPK

image



Background

Alternative Names Suppressor APC domain-containing protein 1 (Protein G7d)

Note

For Research Use Only , Not for Diagnostic Use