


Summary

Name	TIAM2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IVF5
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GKGTLLDSDE GSLSSGTQSS GCPTAEGRQD SKSTSPGKYP HPGLADFADN LIKESDILSD EDDDHRQTVK QGSPTKDIEI QFQRLRISED PDVHPEAEQQ PGPESGEGQK GGEQPKLVRG HFCPIKRKAN STKRDRGTLL KAQIRHQSLD SQSENATIDL NSVLEREFSV

image



Background

Alternative Names	T-lymphoma invasion and metastasis-inducing protein 2 (TIAM-2) (SIF and TIAM1-like exchange factor)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human TIAM2 Protein
Catalog #: PEH9917

