


Product Name: Recombinant Human TBC1D19 Protein
Catalog #: PEH9918



Summary

Name	TBC19_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N5T2
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IIAQIVQKLK GSNLYSQLER QAWASLQRPE IKLESLKEDI KEFFKISGWE KKLQNAVYSE LSVFPLPSHP AAPPEHLKEP LVYMRKAQGS WEKRILKSLN SMCTELSIPL ARKRPVGEQK ELLNKWNEMG TDEPDLSLFR PVYAPKDFLE VLINLRNPNY ENGDSLSEFT HLGLIQVPLK VKDIPELKEC FVELGLNIGQ LGIDDSTQVP

image



Background

Alternative Names TBC1 domain family member 19

Note

Product Name: Recombinant Human TBC1D19 Protein
Catalog #: PEH9918



For Research Use Only , Not for Diagnostic Use