


**Product Name: Recombinant Human TEX9 Protein**  
**Catalog #: PEH9922**



## Summary

|                              |                                                                                                                                                                                                                                                                                                               |
|------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Name</b>                  | TEX9_HUMAN                                                                                                                                                                                                                                                                                                    |
| <b>Purity</b>                | Greater than 85% as determined by reducing SDS-PAGE                                                                                                                                                                                                                                                           |
| <b>Endotoxin level</b>       | Please contact with the lab for this information                                                                                                                                                                                                                                                              |
| <b>Tag</b>                   | C-His                                                                                                                                                                                                                                                                                                         |
| <b>Accession #</b>           | Q8N6V9                                                                                                                                                                                                                                                                                                        |
| <b>Host</b>                  | E.coli                                                                                                                                                                                                                                                                                                        |
| <b>Species</b>               | Human                                                                                                                                                                                                                                                                                                         |
| <b>Molecular Mass</b>        | 33 kDa                                                                                                                                                                                                                                                                                                        |
| <b>Buffer</b>                | Phosphate buffered saline                                                                                                                                                                                                                                                                                     |
| <b>Shipping</b>              | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.                                                                                                                                                                                            |
| <b>Stability&amp;Storage</b> | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.                                                                                                |
| <b>Sequence</b>              | DPVNVQVQNKL HSANKGRKTN SSVKLKYSVD QTADDVAIPE DFSDFSLAKT<br>ISKIEGQLEE EGLPEYIDDI FSGVSNIDGT EAQIRFLKAK LHVMQEELDN<br>VVCECNKKED EIQNLKSQVK NFEEDFMRQQ RTINMQQSQV EKYKTLFEEA<br>NKKYDGLQQQ LSSVERELEN KRRLLQQAAS SQSATEVRLN RALEEAEKYK<br>LELSKLRQNN KDIANEEHKK IEVLKSENKK LEKQKGELMI GFKKQLKLID<br>VLKRQKMHIE |

## image



## Background

|                          |                            |
|--------------------------|----------------------------|
| <b>Alternative Names</b> | Testis-expressed protein 9 |
|--------------------------|----------------------------|

**Product Name: Recombinant Human TEX9 Protein**  
**Catalog #: PEH9922**

---



**Note**

For Research Use Only , Not for Diagnostic Use