


Product Name: Recombinant Human TSSK1B Protein
Catalog #: PEH9923



Summary

Name	TSSK1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BXA7
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SVAINKEGES SRGTEPLWTP EPGSDKKSAT KLEPEGEAQP QAQPETKPEG TAMQMSRQSE ILGFPSKPST METEEGPPQQ

image



Background

Alternative Names	Testis-specific serine/threonine-protein kinase 1 (TSK-1) (TSK1) (TSSK-1) (Testis-specific kinase 1) (EC 2.7.11.1) (Serine/threonine-protein kinase 22A)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use