



Summary

Name	TTC1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q99614
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PEDLLNGLKV TDTQEAECAG PPVPDPKNQH SQSKLLRDDE AHLQEDQGEE ECFHDCSASF EEEPGADKVE NKSNEDEVNSS ELDEEYLIEL EKNMSDEEKQ KRREESTRLK EEGNEQFKKG DYIEAESSYS RALEMCPSCF QKERSILFSN RAAARMKQDK KEMAINDCSK AIQLNPSYIR AILRRAELYE KTDKLDEALE DYKSILEKDP SIHQAREACM RLPKQIEERN ERLKEEMLGK LKDLGNLVLR PFGSLSTENFQ IKQDSSTGSY

image



Background

Alternative Names Tetratricopeptide repeat protein 1 (TPR repeat protein 1)

Product Name: Recombinant Human TTC1 Protein
Catalog #: PEH9925



Note

For Research Use Only , Not for Diagnostic Use