


Product Name: Recombinant Human TTC14 Protein
Catalog #: PEH9927



Summary

Name	TTC14_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96N46
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NKALEIDKQN VEALVARGAL YATKGSLNKA IEDFELALEN CPTHRNARKY LCQTLVERGG QLEEEKFLN AESYKKALA LDETFKDAED ALQKLHKYMQ KSLELREKQA EKEEKQTKK IETSAEKLK LLKEEKRLKK KRRKSTSSSS VSSADESVSS SSSSSSSGHK

image



Background

Alternative Names Tetratricopeptide repeat protein 14 (TPR repeat protein 14)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human TTC14 Protein
Catalog #: PEH9927

