


Product Name: Recombinant Human TTC30B Protein
Catalog #: PEH9928



Summary

Name	TT30B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N4P2
Host	E.coli
Species	Human
Molecular Mass	55 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KACLYAEATR VAFLLLDNPA YHSRVLRLQA AIKYSEGDLP GSRSLVEQLP SREGGEESGG ENETDQQINL GCLLYKEGQY EAACSKFFAA LQASGYQPDL SYNLALAYYS SRQYASALKH IAEIIEGIR QHPGLGVGMT TEGIDVRSVG NTLVLHQTAL VEAFLNKAAL EYQLRNYEAA QEALTMPPR AEEELDPVTL HNQALMNMDA RPTGEFEKLQ FLLQQNPFPP ETFGNLLLLY CKYEYFDLAA DVLAENAHLI YKFLTPYLYD FLDVITCQT APEEAFIKLD GLAGMLTEVL RKLTIQVQEA RHNDRDEAIK KAVNEYDETM EKYPVLMAQ AKIYWNLENY PMVEKIFRKS VEFCDHDVW KLNVAHVLFM QENKYKEAIG FYEPIVKKHY DNILNVAIV LANLCVSYIM TSQNEEAEEL MRKIEKEEQ LSYDDPKKM YHLCIVNLVI GTLYCAKGNV

image



Background

Product Name: Recombinant Human TTC30B Protein
Catalog #: PEH9928



Alternative Names

Tetratricopeptide repeat protein 30B (TPR repeat protein 30B)

Note

For Research Use Only , Not for Diagnostic Use