


Product Name: Recombinant Human TTC8 Protein
Catalog #: PEH9929



Summary

Name	TTC8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TAM2
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LEKSPYDQEP DPELPVHQAA WILKARALTE MVIYIDEIDVD QEGIAEMMLD ENAIQVPRP GTSLKLPGTN QTGGPSQAVR PITQAGRPIT GFLRPSTQSG RPGTMEQAIR TPRTAYTARP ITSSSGRFVR LGTASMLTSP DGPFINLSRL NLTKYSQKPK LAKALFEYIF HHENDVKTIH LEDVVLHLGI YPFLLRNKNH IEKNALDLAA LSTEHSQYKD WWWKVQIGKC YYRLGMYREA EKQFKSALKQ

image



Background

Alternative Names Tetratricopeptide repeat protein 8 (TPR repeat protein 8) (Bardet-Biedl syndrome 8 protein)

Note

Product Name: Recombinant Human TTC8 Protein
Catalog #: PEH9929



For Research Use Only , Not for Diagnostic Use