


Product Name: Recombinant Human TRIP6 Protein
Catalog #: PEH9932



Summary

Name	TRIP6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	Q15654
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GSLDAEIDLL SSTLAELNGG RGHASRRPDR QAYEPPPPPA YRTGSLKPNP ASPLPASPYG GPTPASYTTA STPAGPAFPV QVKVAQPVRG CGPPRRGASQ ASGPLPGPHF PLPGRGEVWG PGYRSQREPG PGAKEEAAGV SGPAGRGRGG EHGPQVPLSQ PPEDELDRLT KKLVDHDMNHP

image



Background

Alternative Names Thyroid receptor-interacting protein 6 (TR-interacting protein 6) (TRIP-6) (Opa-interacting protein 1) (OIP-1) (Zyxin-related protein 1) (ZRP-1)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human TRIP6 Protein
Catalog #: PEH9932

