


Product Name: Recombinant Human TOGARAM2 Protein
Catalog #: PEH9934



Summary

Name	TGRM2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6ZUX3
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KVLVPVAVYC GSIPTSAGP RVLPPGSINS SLPHGEGSLQ PEPRALLNNE EPSQLLRGLG QLGGLKLDTP SKGWQARNGH PRNLRALSLG DQPLVLLPSP ESEANSVARD TIQIKDKLKK RRLSEGLAAS SRASLDPGGG PQGVPLHSTI PRATSQRLLR VPRPMPLIQS IPTTPEASGV

image



Background

Alternative Names TOG array regulator of axonemal microtubules protein 2 (Crescerin-2)

Note

For Research Use Only , Not for Diagnostic Use



Product Name: Recombinant Human TOGARAM2 Protein
Catalog #: PEH9934
