


Product Name: Recombinant Human TMA16 Protein
Catalog #: PEH9937



Summary

Name	TMA16_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96EY4
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GREKKVIHPY SRKAAQITRE AHKQEKKEKL KNEKALRLNL VGEKLQWFQN HLDPQKKRYS KKDACELIER YLNRFSSSELE QIELHNSIRD RQGRRHCSRE TVIKQTMERE RQQFEGYGLE IPDILNASNL KTFREWDFDL KKLPNIMMRK ICANDAIPKT CKRKTIITVD QDLGELELND ESSDSDEEMT

image



Background

Alternative Names Translation machinery-associated protein 16

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human TMA16 Protein
Catalog #: PEH9937

