


Product Name: Recombinant Human TCHHL1 Protein
Catalog #: PEH9939



Summary

Name	TCHL1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5QJ38
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TDQLSLMQLP GKGDSRNELK VQGPSSKEEK GRATEAQNTL LESLDEDNSA SLKIQLETKE PVTSEEEDES PQELAGEGGD QKSPAKKEHN SSVPWSSLEK QMQRDQEPES VERGAVYSSP LYQYLQEKIL QQTNVTQEEH

image



Background

Alternative Names	Trichohyalin-like protein 1 (Basalin) (Protein S100-A17) (S100 calcium-binding protein A17)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use