


Product Name: Recombinant Human TRIM16L Protein
Catalog #: PEH9942



Summary

Name	TR16L_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q309B1
Host	E.coli
Species	Human
Molecular Mass	41 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RKAQANVMLF LEEKEQAALS QANGIKAHLE YRSAEMEKSK QELETMAAIS NTVQFLECYC KFKNTEDITF PSVYIGLKDK LSGIRKVITE STVHLIQLLE NYKKKLQEFK KEEYDIRTQ VSAIVQRKYW TSKPEPSTRE QFLQYVHDIT FDPDTAHKYL RLQEENRKVT NTPWEHPYP DLPSRFLHWR QVLSQQSLYL HRYFFEVEIF GAGTYVGLTC KGIDQKGEER SSCISGNNFS WSLQWNGKEF TAWYSDMETP LKAGPFWRGL VYIDFPGGIL SFYGVEYDSM TLVHKFACKF SEPVYAIFWL SKKENAIRIV DLGEEPEKPA

image



Background

Alternative Names	Tripartite motif-containing protein 16-like protein (Tripartite motif-containing protein 70)
--------------------------	--

Product Name: Recombinant Human TRIM16L Protein
Catalog #: PEH9942



Note

For Research Use Only , Not for Diagnostic Use