


Product Name: Recombinant Human TRIM73 Protein
Catalog #: PEH9944



Summary

Name	TRI73_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q86UV7
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DGSSSLPNVS LAWVIEALRL PGDPEPKVCV HHRNPLSLFC EKDQELICGL CGLLGSHQHH PVTPVSTVCS RMKEELAALF SELKQEQKKV DELIAKLVKN RTRIVNESDV FSWVIRREFQ ELRHPVDEEK ARCLEGIGGH TRGLVASLDM QLEQAQGTRE RLAQAECVLE

image



Background

Alternative Names Tripartite motif-containing protein 73 (Tripartite motif-containing protein 50B)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human TRIM73 Protein
Catalog #: PEH9944

