


Product Name: Recombinant Human TRMT10C Protein
Catalog #: PEH9946



Summary

Name	TM10C_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q7L0Y3
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SSKIPAVTYP KNESTPPSEE LELDKWKTTM KSSVQEECVS TISSSKDEDP LAATREFIEM WRLLGREVPE HITEEELKTL MECVSNTAKK KYLKLYLTKE KVKKARQIKK EMKAAAREEA KNIKLETTE

image



Background

Alternative Names	tRNA methyltransferase 10 homolog C (HBV pre-S2 trans-regulated protein 2) (Mitochondrial ribonuclease P protein 1) (Mitochondrial RNase P protein 1) (RNA (guanine-9-)-methyltransferase domain-containing protein 1) (Renal carcinoma antigen NY-REN-49) (mRNA methyladenosine-N(1)-methyltransferase) (EC 2.1.1.-) (tRNA (adenine(9)-N(1))-methyltransferase) (EC 2.1.1.218) (tRNA (guanine(9)-N(1))-methyltransferase) (EC 2.1.1.221)
--------------------------	---

Product Name: Recombinant Human TRMT10C Protein
Catalog #: PEH9946



Note

For Research Use Only , Not for Diagnostic Use