


Product Name: Recombinant Human TYW5 Protein
Catalog #: PEH9947



Summary

Name	TYW5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	A2RUC4
Host	E.coli
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LEGVSREQFM QHLYPQRKPL VLEGIDLGPC TSKWTVDYLS QVGGKKEVKI HVAAVAQMDF ISKNFVYRTL PFDQLVQRAA EEKHKEFFVS EDEKYYLRSL GEDPRKDVAD IRKQFPLLKG DIKFPEFFKE EQFFSSVFRI SSPGLQLWTH YDVMDNLLIQ VTGKKRVVLF SPRDAQYLYL KGTKSEVLNI DNPDLAKYPL FSKARRYECS LEAGDVLFI ALWFHNVISE EFGVGVNIFW KHLPSECYDK TDTYGNKDPT AASRAAQILD RALKTLAELP EYRDFYARR MVLHIQDKAY

image



Background

Alternative Names	tRNA wybutosine-synthesizing protein 5 (hTYW5) (EC 1.14.11.42) (tRNA(Phe) (7-(3-amino-3-carboxypropyl)wyosine(37)-C(2))-hydroxylase)
--------------------------	--

Product Name: Recombinant Human TYW5 Protein
Catalog #: PEH9947



Note

For Research Use Only , Not for Diagnostic Use