


Product Name: Recombinant Human DUS2 Protein
Catalog #: PEH9948



Summary

Name	DUS2L_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NX74
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	STVDFVAPDD RLVFRTTCERE QNRVVFQMGT SDAERALAVA RLVENDVAGI DVNMGCPKQY STKGGMGAAL LSDPKIEKI

image



Background

Alternative Names	tRNA-dihydrouridine(20) synthase [NAD(P)+]-like (EC 1.3.1.91) (Dihydrouridine synthase 2) (Up-regulated in lung cancer protein 8) (URLC8) (tRNA-dihydrouridine synthase 2-like) (hDUS2)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use