


**Product Name: Recombinant Human TUBB2B Protein**  
**Catalog #: PEH9951**



## Summary

|                              |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Name</b>                  | TBB2B_HUMAN                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Purity</b>                | Greater than 85% as determined by reducing SDS-PAGE                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Endotoxin level</b>       | Please contact with the lab for this information                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Tag</b>                   | C-His                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Accession #</b>           | Q9BVA1                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Host</b>                  | E.coli                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Species</b>               | Human                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Molecular Mass</b>        | 51 kDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Buffer</b>                | Phosphate buffered saline                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Shipping</b>              | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Stability&amp;Storage</b> | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.                                                                                                                                                                                                                                                                                                    |
| <b>Sequence</b>              | QCGNQIGAKF WEVISDEHGI DPTGSYHGDS DLQLERINVY YNEATGNKYV<br>PRAILVDLEP GTMDSVRS GP FGQIFRPDNF VFGQSGAGNN WAKGHYTEGA<br>ELVDSVLDVV RKESESCDCL QGFQLTHSLG GGTGSGMGTL LISKIREEYP<br>DRIMNTFSVM PSPKVS DTVV EPYNATLSVH QLVENTDETY CIDNEALYDI<br>CFRTLKLTTP TYGDLNHLVS ATMSGVTTCL RFPGQLNADL RKLAVNMVFP<br>PRLHFFMPGF APLTSRGSQQ YRALTVPELT QQMFDSKNMM AACDPRHGRY<br>LTVAAIFRGR MSMKEVDEQM LNVQNKNSY FVEWIPNNVK TAVCDIPPRG<br>LKMSATFIGN STAIQELFKR ISEQFTAMFR RKAFLHWYTG EGMDMEFETE<br>AESNMNDLVS EYQQYQDATA DEQGEFEEEE |

## image



## Background

**Product Name: Recombinant Human TUBB2B Protein**  
**Catalog #: PEH9951**



---

**Alternative Names**

Tubulin beta-2B chain

**Note**

For Research Use Only , Not for Diagnostic Use