


Product Name: Recombinant Human TUBB6 Protein
Catalog #: PEH9952



Summary

Name	TBB6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BUF5
Host	E.coli
Species	Human
Molecular Mass	51 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QCGNQIGTKF WEVISDEHGI DPAGGYVGDS ALQLERINVY YNESSSQKYV PRAALVDLEP GTMDSVRS GP FGQLFRPDNF IFGQTGAGNN WAKGHYTEGA ELVDAVL DVV RKECEHCDCL QGFQLTHSLG GGTGSGMGTL LISKIREEFP DRIMNTFSVM PSPKVS DTVV EPYNATLSVH QLVENTDETY CIDNEALYDI CFRTLKLTTP TYGDLNHLVS ATMSGVTTSL RFPGQLNADL RKLAVNMVFP PRLHFFMPGF APLTSRGSQQ YRALTVP ELT QQMFDARNMM AACDPRHG RY LTVATVFRGP MSMKEVDEQM LAIQSKNSSY FVEWIPNNVK VAVCDIPPRG LKMASTFIGN STAIQELFKR ISEQFSAMFR RKAFLHWFTG EGMDMEFE TE AESNMNDLVS EYQQYQDATA NDGEEAFEDE

image



Background

Product Name: Recombinant Human TUBB6 Protein
Catalog #: PEH9952



Alternative Names

Tubulin beta-6 chain (Tubulin beta class V)

Note

For Research Use Only , Not for Diagnostic Use