


Product Name: Recombinant Human TUBB8 Protein
Catalog #: PEH9953



Summary

Name	TBB8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q3ZCM7
Host	E.coli
Species	Human
Molecular Mass	49 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	WEVISDEHAI DSAGTYHGDS HLQLERINVY YNEASGGRYV PRAVLVDLEP GTMDSVRSGP FGQVFRPDNF IFGQCGAGNN WAKGHYTEGA ELMESVMDVV RKEAESCDC L QGFQLTHSLG GGTGSGMGTL LLSKIREEYP DRIINTFSIL PSPKVSDTVV EPYNATLSVH QLIENADETF CIDNEALYDI CSKTLKLPTP TYGDLNHLVS ATMSGVTTCL RFPGQLNADL RKLAVNMVPF PRLHFFMPGF APLTSRGSQQ YRALTVAELT QQMFDAKNMM AACDPRHG RY LTAAAI FRGR MPMREVDEQM FNIQDKNSSY FADWLPNNVK TAVCDIPPRG LKMSATFIGN NTAIQELFKR VSEQFTAMFR RKAFLHWYTG EGMDMEFTE AESNMNDLVS EYQQYQDATA EEEEEDEYAE

image



Background

Product Name: Recombinant Human TUBB8 Protein
Catalog #: PEH9953



Alternative Names

Tubulin beta-8 chain (Tubulin beta 8 class VIII)

Note

For Research Use Only , Not for Diagnostic Use