


Product Name: Recombinant Human TTL13 Protein
Catalog #: PEH9955



Summary

Name	TTL13_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	A6NNM8
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ETAASKAREE CARQLEEIR LKQEQQETSG TKRQKARDQN QGESAGEKSR PRAGLQSLST HLAYRNRNWE KELLPGQLDT MRPQEIVEEE ELERMKALLQ RETLIRSLGI

image



Background

Alternative Names Tubulin polyglutamylase TTL13P (EC 6.-.-.-) (Tubulin tyrosine ligase-like family member 13 pseudogene) (Tubulin--tyrosine ligase-like protein 13)

Note

For Research Use Only , Not for Diagnostic Use