


Product Name: Recombinant Human USP32 Protein
Catalog #: PEH9960



Summary

Name	UBP32_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NFA0
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KRYWLLKAQS RTGRFDLETG GPLVSPPIRP SLSEGLFNAF DENRDNHIDF KEISCGLSAC CRGPLAERQK FCFKVFVDVDR DGVLSRVELR DMVVALLEVW KDNRTDDIPE LHMDLSDIVE GILNAHDTTK MGHILTLEDYQ

image



Background

Alternative Names	Ubiquitin carboxyl-terminal hydrolase 32 (EC 3.4.19.12) (Deubiquitinating enzyme 32) (Renal carcinoma antigen NY-REN-60) (Ubiquitin thioesterase 32) (Ubiquitin-specific-processing protease 32)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human USP32 Protein
Catalog #: PEH9960

