


Product Name: Recombinant Human USP33 Protein
Catalog #: PEH9961



Summary

Name	UBP33_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TEY7
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EKRRKTELEI FIRLNRAFQK EDSPATFYCI SMQWFREWES FVKGKGDGPP GPIDNTKIAV TKCGNVMLRQ GADSGQISEE TWNFLQSIYG GGPEVILRPP VVHVDPDILQ AEEKIEVETR

image



Background

Alternative Names Ubiquitin carboxyl-terminal hydrolase 33 (EC 3.4.19.12) (Deubiquitinating enzyme 33) (Ubiquitin thioesterase 33) (Ubiquitin-specific-processing protease 33) (VHL-interacting deubiquitinating enzyme 1) (hVDU1)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human USP33 Protein
Catalog #: PEH9961

