


Summary

Name	UBP42_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H9J4
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IKNPPHLNGT GPLKDTPESS MSSPNGNSSV NRASPVNASA SVQNWSVNRS SVIPEHPKKQ KITISIHNLK PVRQCQSQPN LHSNSLENPT KPVPSTITN SAVQSTSNAS TMSVSSKVTK PIPRSESCSQ PVMNGKSKLN SSVLVPGAE SSEDSDEESK GLGKENGIGT IVSSHSPGQD AEDEEATPHE LQEPMTLNGA NSADSDSDPK ENGLAPDGAS CQGQPALHSE NPFKANGLP GKLMPPALLS LPEDKILETF RLSNKLKGST DEMSAPGAER

image



Background

Alternative Names	Ubiquitin carboxyl-terminal hydrolase 42 (EC 3.4.19.12) (Deubiquitinating enzyme 42) (Ubiquitin thioesterase 42) (Ubiquitin-specific-processing protease 42)
--------------------------	--

Product Name: Recombinant Human USP42 Protein
Catalog #: PEH9963



Note

For Research Use Only , Not for Diagnostic Use