


Product Name: Recombinant Human UBFD1 Protein
Catalog #: PEH9965



Summary

Name	UBFD1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O14562
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RGSLQPAPAQ PPGDPAAQAS VSNGEDAGGG AGRELVDLKI IWNKTKHDVK FPLDSTGSEL KQKIHSITGL PPAMQKVMYK GLVPEDKTLR EIKVTSGAKI MVGSTINDV LAVNTPKDAA QQDAKAEENK KEPLCRQKQH RKVLDKGKPE DVMPSVKGAQ

image



Background

Alternative Names Ubiquitin domain-containing protein UBFD1 (Ubiquitin-binding protein homolog)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human UBFD1 Protein
Catalog #: PEH9965

