


Product Name: Recombinant Human C10orf67 Protein
Catalog #: PEH9970



Summary

Name	CJ067_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IYJ2
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FFSTDHATQT DSSEILSVKE LSSSTQKLAQ MMKSLQVDFG FLKQLLQLKF EDRLKEESLS LFTILHDRIL EIEKHYQQNE DKMRKSFNQQ LADAIKVIK MYQQFFEVEE ENVSLQDAST VKTNILLRKL KEKEEVIKEL KEELDQYKDF GFHKMESFAK ETSSPKSNLE KENLEYKVEN ERLQLISEL EEEIQINLKE NSGLEDELIS

image



Background

Alternative Names	Uncharacterized protein C10orf67, mitochondrial
--------------------------	---

Note

Product Name: Recombinant Human C10orf67 Protein
Catalog #: PEH9970



For Research Use Only , Not for Diagnostic Use