


Product Name: Recombinant Human C1orf141 Protein
Catalog #: PEH9976



Summary

Name	CA141_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5J VX7
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AISKIKEDKS CSITKSKMHV SFKCEPEPRK SNFEKSNLRP FFIQTNVKNK ESESTAQIEK KPRKPLDSVG LLEGDRNKRK KSPQMND FNI KENKSVRNYQ LSKYRSVRKK SLLPLCFEDE LKNPHAKIVN VSPTKTVTSH MEQKDTNP II

image



Background

Alternative Names Uncharacterized protein C1orf141

Note

For Research Use Only , Not for Diagnostic Use