


Product Name: Recombinant Human C4orf19 Protein
Catalog #: PEH9977



Summary

Name	CD019_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IY42
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IGPAATPQPT GNSPTQDDR GSWASTANTV PPTQPFLEGG GTRKQDCVLL ASEGTQVMRN GDSRAPSEAE SFALEVQDHV FQIPAPDYLQ HWGPAGDNVD HNEKDCVFKN HTEDESLEGI QPPVGEHGLN TPFSVRRSWD SLNEDVETEV LSICFNEKGP VHAMPVVDSG NRQEDTHGSD GDGDGEIVDE DAAVAEALAA LEAATAGEDL

image



Background

Alternative Names	Uncharacterized protein C4orf19
--------------------------	---------------------------------

Note

Product Name: Recombinant Human C4orf19 Protein
Catalog #: PEH9977



For Research Use Only , Not for Diagnostic Use