


Product Name: Recombinant Human C7orf50 Protein
Catalog #: PEH9978



Summary

Name	CG050_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BRJ6
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VTEKKNKKLK KASAEGPLLG PEAAPSGEGA GSKGEAVLRP GLDAEPELSP EEQRVLERKL KKERKKEERQ RLREAGLVAQ HPPARRSGAE LALDYLCRWA QKHKNWRFQK TRQTWLLHM YSDKVPDEH FSTLLAYLEG LQGRARELTV QKAEALMREL DEEGSDPPLP

image



Background

Alternative Names Uncharacterized protein C7orf50

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human C7orf50 Protein
Catalog #: PEH9978

