



Product Name: Recombinant Human ZSWIM9 Protein
Catalog #: PEH9980



Summary

Name	ZSWIM9_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q86XI8
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GAQGENEVR GLETGDWGGG PKEGSIWRGA QMEKEWARAL ETRDWGGAQF EGEKGRALQI RDWRGGRLN QKPRGLEGGV LRGSKLEKGH LRGPEIRDWR GPQLEGEKDW GLEGYVWRGS QLEDQALRGL EGYTWRVAQL EDRRSTDLR GTQFDYER GK

image



Background

Alternative Names Uncharacterized protein ZSWIM9

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ZSWIM9 Protein
Catalog #: PEH9980

