


Product Name: Recombinant Human PHAF1 Protein
Catalog #: PEH9983



Summary

Name	CP070_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BSU1
Host	E.coli
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SLGNEQWEFT LGMPLAQAVA ILQKHCRIIK NVQVLYSEQS PLSHDLILNL TQDGIKLMFD AFNQRLKVIE VCDLTKVKLK YCGVHFNSQA IAPTIEQIDQ SFGATHPGVY NSAEQLFHLN FRGLSFSFQL DSWTEAPKYE PNFAHGLASL QIPHGATVKR MYIYSGNSLQ DTKAPMMPLS CFLGNVYAES VDVLRDGTGP AGLRLRLAA GCGPGLLADA KMRVFERSVY FGDSCQDVLS MLGSPHKVVFY KSEDKMKIHS PSPHKQVPSK CNDYFFNYFT

image



Background

Alternative Names UPF0183 protein C16orf70

Product Name: Recombinant Human PHAF1 Protein
Catalog #: PEH9983



Note

For Research Use Only , Not for Diagnostic Use