


Product Name: Recombinant Human C5orf22 Protein
Catalog #: PEH9986



Summary

Name	CE022_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q49AR2
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TAASTGEILE ILKKGKAFVL DIDLDFFSVK NPFKEMFTQE EYKILQELYQ FKKPGTNLTE EDLVDIVDTR IHQLEDLEAT FADLCDGDDE ETVQRWASNP GMESLVPLVQ SLKKRMEVPD YEMVHQAGLT CDYSELPHHI STEQEIECLI QSVHYLLKNL PNPTLVTIAR SSLDDYCPSD QVDTIQEKVL

image



Background

Alternative Names UPF0489 protein C5orf22

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human C5orf22 Protein
Catalog #: PEH9986

