


Product Name: Recombinant Human VPS11 Protein
Catalog #: PEH9989



Summary

Name	VPS11_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q9H270
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KQSQQIAQDE LRVRRYREET TRIRQEIQL KASPKIFQKT KCSICNSALE LPSVHFLCGH SFHQHCFESY SESDADCPTC LPENRKVMDM IRAQEQKRDL HDQFQHQLRC SNDSFSVIAD YFGRGVFNKL TLLDPPTAR LTSSLEAGLQ RDLLMHSRRG

image



Background

Alternative Names	Vacuolar protein sorting-associated protein 11 homolog (hVPS11) (RING finger protein 108)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human VPS11 Protein
Catalog #: PEH9989

