


Product Name: Recombinant Human VPS39 Protein
Catalog #: PEH9991



Summary

Name	VPS39_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q96JC1
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DKQFELALQL AEMKDDSDSE KQQIHHIKN LYAFNLCQK RFDESMQVFA KLGTDPTHVM GLYPDLLPTD YRKQLQYPNP

image



Background

Alternative Names Vam6/Vps39-like protein (TRAP1-like protein) (hVam6p)

Note

For Research Use Only , Not for Diagnostic Use