



Product Name: Recombinant Human WDR7 Protein
Catalog #: PEH9996



Summary

Name	WDR7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y4E6
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LFFDVEALII QLLTEEASRP NTALISPENL QKASGSSDKG GSFLTGKRAA VLFQQVKETI KENIKEHLLD DEEEDDEEIMR QRREESDPEY RSSKSKPLTL LEYNLTMDDTA KLFMSCLHAW GLNEVLDEVC LDRLGMLKPH CTVSFGLLSR GGHMSLMLPG YNQPACKLSH GKTEVGRKLP

image



Background

Alternative Names WD repeat-containing protein 7 (Rabconnectin-3 beta) (TGF-beta resistance-associated protein TRAG)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human WDR7 Protein
Catalog #: PEH9996

