


Product Name: Recombinant Human YBX3 Protein
Catalog #: PEH9998



Summary

Name	YBOX3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P16989
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	APAATGTAAA ASLATAAGSE DAEKKVLATK VLGTVKWFNV RNYGYGFNIRN DTKEDVFVHQ TAIKKNNPRK YLRVGDGET VEFDVVEGEK GAEEAANVTGP DGVPEGSRY AADRRRYRRG YYGRRRGPPR NYAGEEEEEEG SGSSEGFDP ATDRQFSGAR NQLRRPQYRP QYRQRRFPY HVGQTFDRRS RVLPHPNRIQ AGEIGEMKDG

image



Background

Alternative Names	Y-box-binding protein 3 (Cold shock domain-containing protein A) (DNA-binding protein A) (Single-strand DNA-binding protein NF-GMB)
--------------------------	---

Note

Product Name: Recombinant Human YBX3 Protein
Catalog #: PEH9998



For Research Use Only , Not for Diagnostic Use