


Product Name: Recombinant Mouse Beta-NGF (110AA)
Catalog #: PEM0151



Summary

Name	beta-NGF/ β -Nerve Growth Factor/ β -NGF (Met130-Arg239)
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/ μ g as determined by LAL test.
Construction	Recombinant Mouse Beta-Nerve Growth Factor is produced by our E.coli expression system and the target gene encoding Met130-Arg239 is expressed.
Accession #	P01139
Host	E.coli
Species	Mouse
Predicted Molecular Mass	12.4 KDa
Formulation	Lyophilized from a 0.2 μ m filtered solution of 20mM PB, 200mM NaCl, pH 8.0.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Lyophilized protein should be stored at $\leq -20^{\circ}\text{C}$, stable for one year after receipt. Reconstituted protein solution can be stored at 2-8 $^{\circ}\text{C}$ for 2-7 days. Aliquots of reconstituted samples are stable at $\leq -20^{\circ}\text{C}$ for 3 months.
Reconstitution	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μ g/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μ g/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE image



Background

Product Name: Recombinant Mouse Beta-NGF (110AA)
Catalog #: PEM0151



Alternative Names

Beta-Nerve Growth Factor; Beta-NGF; NGF; NGFB; β -NGF

Background

NGF is the first member discovered in the Neurotrophin family, which includes brain-derived neurotrophic factor (BDNF), neurotrophin-3 (NT-3), and neurotrophin-4 (NT-4). These proteins belong to the cysteine-knot family of growth factors that assume stable dimeric structures. Mouse beta -NGF is a homodimer of two 120 amino acid polypeptides. It shares approximately 90% homology at the amino acid level with human beta -NGF and 95.8% with rat beta -NGF. NGF signaling has been shown to play an important role in neuroprotection and repair. β -NGF acts as a growth and differentiation factor for B lymphocytes, and enhances B-cell survival. It is a potent neurotrophic factor that signals through its receptor β -NGFR, and plays a crucial role in the development and preservation of the sensory and sympathetic nervous systems.

Note

For Research Use Only , Not for Diagnostic Use.