


Product Name: Recombinant Mouse Galectin-9 (N-GST)
Catalog #: PEM0708



Summary

Name	Galectin-9/LGALS9
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Mouse Galectin-9 is produced by our E.coli expression system and the target gene encoding Met1-Thr322 is expressed with a GST tag at the N-terminus.
Accession #	O08573-2
Host	E.coli
Species	Mouse
Predicted Molecular Mass	62.8 KDa
Formulation	Supplied as a 0.2 μm filtered solution of 20mM PB, 1mM EDTA, 5mM β-Mercaptoethanol, 5% Trehalose, 400mM NaCl, pH6.5.
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Reconstitution	

SDS-PAGE image



Background

Alternative Names	Galectin-9; Lgals9; Gal-9
Background	Galectin-9 is a cytoplasmic protein that contains two galectin domains. Galectin-9

Product Name: Recombinant Mouse Galectin-9 (N-GST)
Catalog #: PEM0708

is an S-type lectin that is over-expressed in Hodgkins disease tissue. Galectin-9 binds galactosides and has high affinity for the Forssman pentasaccharide. Galectin-9 plays a role in thymocyte-epithelial interactions relevant to the biology of the thymus and Inhibits cell proliferation. Galectin-9 is a ligand for HAVCR2/TIM3 and induces T-helper type 1 lymphocyte (Th1) death. In addition, Galectin-9 suppresses tumor cell metastasis by interfering with the associations CD44, VCAM-1, Integrin $\alpha 4\beta 1$.

Note

For Research Use Only , Not for Diagnostic Use.