

Product Name: Recombinant Mouse S100A4 (C-6His)
Catalog #: PEM1456




Summary

Name	S100A4
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Mouse Protein S100-A4 is produced by our E.coli expression system and the target gene encoding Met1-Lys101 is expressed with a 6His tag at the C-terminus.
Accession #	P07091
Host	E.coli
Species	Mouse
Predicted Molecular Mass	12.5 KDa
Formulation	Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Lyophilized protein should be stored at ≤ -20°C, stable for one year after receipt. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at ≤ -20°C for 3 months.
Reconstitution	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE image

Product Name: Recombinant Mouse S100A4 (C-6His)
Catalog #: PEM1456



Alternative Names

Protein S100-A4; Metastasin; Metastatic cell protein; PEL98; Placental calcium-binding protein; Protein 18A2; Protein Mts1; S100 calcium-binding protein A4; S100a4; Capl; Mts1

Background

S100A4 is a member of the S100 family of proteins. The S100 family is further classified as a member of the EF-hand superfamily of Ca⁺⁺-binding proteins. These participate in both calcium-dependent and calcium-independent protein-protein interactions. The hallmark of this superfamily is the EF-hand motif that consists of a Ca⁺⁺-binding site flanked by two α -helices (helix E and helix F) that were originally identified in a right-handed model of carp muscle calcium-binding protein. Mouse S100A4 is 101 amino acids (aa) in length. It contains two EF hand domains, one between aa 12-47, and a second between aa 50-85. S100A4 activity has been associated with cell transformation. It seems likely this is either coincidental, or a consequence, rather than a cause of transformation.

Note

For Research Use Only , Not for Diagnostic Use.