


Summary

Name	Alpha 1-Microglobulin/AMBP
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Human Alpha 1-Microglobulin is produced by our Mammalian expression system and the target gene encoding Gly20-Val203 is expressed with a 6His tag at the C-terminus.
Accession #	P02760
Host	Human Cells
Species	Human
Predicted Molecular Mass	21.9 KDa
Formulation	Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Lyophilized protein should be stored at ≤ -20°C, stable for one year after receipt. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at ≤ -20°C for 3 months.
Reconstitution	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE image

Product Name: Recombinant Human AMBP (C-6His)
Catalog #: PHH0045



Alternative Names

Protein AMBP; Alpha-1-Microglobulin; Protein HC; Alpha-1 Microglycoprotein; Complex-Forming Glycoprotein Heterogeneous in Charge; Inter-Alpha-Trypsin Inhibitor Light Chain; ITI-LC; Bikunin; EDC1; HI-30; Uronic-Acid-Rich Protein; Trypstatin; AMBP; HCP; ITIL

Background

Protein AMBP belongs to the calycin superfamily and Lipocalin family. AMBP can be cleaved into three chains: α -1-microglobulin, inter- α -trypsin inhibitor light chain and trypstatin. AMBP is expressed by the liver and secreted in plasma. α -1-microglobulin occurs in many physiological fluids including the plasma, urine, and cerebrospinal fluid. Inter- α -trypsin inhibitor is present in the plasma and urine. α -1-microglobulin occurs as a monomer and also in complexes with IgA and albumin, Inter- α -trypsin inhibitor inhibits trypsin, plasmin and lysosomal granulocytic elastase. Trypstatin act as a trypsin inhibitor, exists in a monomer forms and also occurs as a complex with trypstatin in mast cells.

Note

For Research Use Only , Not for Diagnostic Use.